

A proactive breakthrough in cancer screening



Galleri[®] is a proactive blood test that can be done annually, screening for many deadly cancers before symptoms appear. With a single blood draw, the Galleri test screens for a “fingerprint” of 50+ cancers, including many without recommended screening available today (like pancreatic, ovarian, liver, and more).*

The Galleri test does not detect a signal for all cancers and not all cancers can be detected in the blood. False positive and false negative results do occur. Galleri is a screening test and does not diagnose cancer. Diagnostic testing is needed to confirm cancer. The Galleri test should be used in addition to healthcare provider recommended screening tests. Galleri is available by prescription only.

*Sensitivity in study participants — Pancreatic cancer: 83.7% overall (61.9% stage I, 60.0% stage II, 85.7% stage III, 95.9% stage IV); Ovarian cancer: 83.1% overall (50.0% stage I, 80.0% stage II, 87.1% stage III, 94.7% stage IV); Liver & Bile-ducts cancer: 93.5% overall (100% stage I, 70.0% stage II, 100% stage III, 100% stage IV).

Galleri for Granite

The Galleri test is available to employees and dependents enrolled in an Anthem medical plan ages 50 years or older, or 40-49 with elevated risk factors.

The Galleri test is not recommended if you are pregnant or undergoing active cancer treatment.

While individuals ages 40-49 with certain risk factors may be at an elevated risk for cancer, Galleri was only studied in a population ages 50 and older.



Important Safety Information The Galleri test is recommended for use in adults with an elevated risk for cancer, such as those aged 50 or older. The Galleri test does not detect all cancers and should be used in addition to routine cancer screening tests recommended by a healthcare provider. Galleri is intended to detect cancer signals and predict where in the body the cancer signal is located. Use of Galleri is not recommended in individuals who are pregnant, 21 years old or younger, or undergoing active cancer treatment. Results should be interpreted by a healthcare provider in the context of medical history, clinical signs and symptoms. A test result of No Cancer Signal Detected does not rule out cancer. A test result of Cancer Signal Detected requires confirmatory diagnostic evaluation by medically established procedures (e.g. imaging) to confirm cancer. If cancer is not confirmed with further testing, it could mean that cancer is not present or testing was insufficient to detect cancer, including due to the cancer being located in a different part of the body. False-positive (a cancer signal detected when cancer is not present) and false-negative (a cancer signal not detected when cancer is present) test results do occur. **Rx only.**

Laboratory/Test Information The GRAIL clinical laboratory is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and accredited by the College of American Pathologists. The Galleri test was developed and its performance characteristics were determined by GRAIL. The Galleri test has not been cleared or approved by the Food and Drug Administration. The GRAIL clinical laboratory is regulated under CLIA to perform high-complexity testing. The Galleri test is intended for clinical purposes.



Learn more about the Galleri test and the benefit at Galleri.com/Granite

Simple process



Request the test

It will take about 10 minutes to complete a health questionnaire. Galleri is available by prescription only. By requesting the test and completing the questionnaire, you're requesting a prescription from an independent telemedicine provider.



Receive a collection kit

If the telemedicine provider confirms your eligibility and prescribes Galleri, you will receive a collection kit in the mail.



Provide your blood sample

A trained technician will draw your blood sample and ship the kit to the GRAIL lab.



Receive your results

You'll get your results about 2 weeks after your blood is received at the GRAIL lab.

How does Galleri work?

Cancers growing in the body shed DNA into the bloodstream. Although there are many types of cancer, the DNA fragments act like a unique "fingerprint" of cancer. With this unique "fingerprint" of cancer, the Galleri test helps to provide direction to a healthcare provider on the cancer's origin and the next steps in diagnosis.

Understanding your results

Your result will give you the confidence to know you've done more to screen for cancer. There are two possible test results:

**! Cancer Signal Detected +
Cancer Signal Origin prediction**

EVALUATION FOR CANCER

Diagnostic evaluation is needed to confirm cancer

✓ No Cancer Signal Detected

ROUTINE CARE

Continue recommended screening + consider annual Galleri testing

The Galleri test does not detect a signal for all cancers and not all cancers can be detected in the blood. False positive and false negative results do occur. The Galleri test identifies DNA in the bloodstream shed by cancer cells and does not predict future genetic risk for cancer. Galleri should be used in addition to healthcare provider recommended screening tests.



Check eligibility and request the Galleri test at [Galleri.com/Granite](https://www.Galleri.com/Granite)

Galleri is available by prescription only. By requesting the test and completing the health questionnaire, you're requesting a prescription from an independent healthcare provider. The provider will review your questionnaire and prescribe Galleri if appropriate.